

Upper Arm

Digital Blood Pressure Monitor

Take care of your health



One click measurement



HD screen with backlight



2*99 groups of memory



self-measure with cuff



Voice guidance



Remind arrhythmia



Display blood pressure



Divide into A/B groups

One-click
Measurement

Voice
Announcements

Irregular
Pulse

90 Groups
of memory

Dual power
supply mode

1

One click start Measure automatically

One key measurement, easy to operate
and then the cuff will tighten up and the machine will measure automatically



One-click
Measurement

Voice
Announcements

Irregular
Pulse

90 Groups
of memory

Dual power
supply mode

Voice prompt Auto broadcast

The whole process voice prompts the user how to use the machine, and automatically broadcast the results after measurement, which is convenient for the elderly



Your systolic blood
pressure is 120mmhg
.....

Please keep
silence and relax



One-click
Measurement

Voice
Announcements

Irregular
Pulse

90 Groups
of memory

Dual power
supply mode

3

Heart rate monitoring

WHO Hypertension warning

With irregular heart rate prompt, this  shape will appear when detecting irregular heart rate, timely monitoring irregular heart rate, and voice broadcast reminder when detecting hypertension

180
mmHg

100
mmHg


Abnormal

80/min

Irregular heart rate
symbol



One-click
Measurement

Voice
Announcements


Irregular
Pulse

90 Groups
of memory

Dual power
supply mode

2*99 groups of memory

It can record 90 groups of measurement data of two users, easily help parents record the blood pressure value, so that children can observe their parents' health at any time



One-click
Measurement

Voice
Announcements

Irregular
Pulse

90 Groups
of memory

Dual power
supply mode


5

Dual power supply mode with batteries and power cord

There are two power supply modes that it can powered with batteries as well as charged with powered cord. So it is convenient for the users to charge the machine.



How to use it correctly



Wear the cuff

According to the figure, take off the thicker clothes of the measured arm, Arms up, flat on the table, Wear the cuff to 2-3cm of the upper arm.




Sitting posture

Sit in a chair and relax with the center of the cuff at the same height as the heart, wait for a minutes and it will come out the result immediately.



Panel display



Product Display

